


# Certified Reference Material (CRM)

- Certified Reference Powders dedicated to Particles Size Analysers **Calibration & Verification**
- Specific kit for **Dry & Liquid** dispersion mode
- Essential for Quality System in the frame of the **ISO13320** norm
- **Traceable** to international standards

- Delivered in **Single Shot vials** packaged in a specific box

- Complete set of **certificates and procedures**




# Technical specification

## Cilas CSR

Cilas **S**econdary **R**eference Powders are used for laboratory calibration:

- Quartz Powders retraced to international standards
- Single shot vials : 0.04gr, 0.233gr, 0.12gr, 0.7gr
- Delivered with **Traceability Certificate**



## WhiteHouse Scientific CRM

White House materials are used for ISO **verification** of the instrument

- Certified glass beads
- Single shot vials : 0.5gr, 1.0gr

## Ordering Information

Verification pack for the 930 or 1064 liquid model	ref: 1061018
Verification pack for the 1180 liquid model	ref: 1061020
Verification pack for the 930, 1064 or 1180 dry model device	ref: 1061030
Calibration pack for the 930 or 1064 liquid model instruments	ref: 1061017
Calibration pack for the 1180 liquid model instrument	ref: 1061019
Calibration pack for the 930, 1064 or 1180 dry dispersion unit	ref: 1061029

More info: [www.particle-size-analyzer.com](http://www.particle-size-analyzer.com)

[info.particle@cilas.com](mailto:info.particle@cilas.com)

Your partner in particle analysis